

Anti-Cox Iv Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P13073

Item Number: CM0003614



Introduction

High-quality rabbit monoclonal antibody against human COX IV (cytochrome c oxidase subunit 4 isoform 1). Validated for WB, IHC-P, and ELISA; cross-reacts with human, mouse, and rat. Ideal for mitochondrial research and diagnostics.

[Learn More](#)

Parameter	Value
Product Name	COX IV Rabbit mAb
Remarks/Alias	COX4; COXIV; COX4-1; COXIV-1; MC4DN16; COX IV-1; COX IV
Species	Human
Gene ID (Human)	1327
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13073
Protein Weight	20kDa
Shipping	Ice bag