


## Anti-Cox Iv Rabbit Monoclonal Antibody For Wb, Ip, If/Icc, Ihc-P, Elisa - P13073

Item Number: CM0004668



### Introduction

Rabbit monoclonal antibody against COX IV (COX4I1), validated for western blot, IP, immunofluorescence, IHC, and ELISA. Reacts with human, mouse, rat. Ideal for mitochondrial respiration and oxidative phosphorylation studies.

[Learn More](#)

| Parameter        | Value                                                                                                                                  |
|------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | [KD Validated] COX IV Rabbit mAb                                                                                                       |
| Remarks/Alias    | COX4; COXIV; COX4-1; COXIV-1; MC4DN16; COX IV-1                                                                                        |
| Species          | Human                                                                                                                                  |
| Gene ID (Human)  | 1327                                                                                                                                   |
| Gene ID          | 1327                                                                                                                                   |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 23-98 of human COX IV. (NP_001852.1) |
| Source           | Rabbit                                                                                                                                 |
| Category         | Monoclonal Antibodies                                                                                                                  |
| Application      | WB, IP, IF/ICC, IHC-P, ELISA                                                                                                           |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                      |
| SWISS            | P13073                                                                                                                                 |
| Protein Weight   | 20 kDa                                                                                                                                 |
| Shipping         | Ice bag                                                                                                                                |