


Anti-Cox1 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P00397

Item Number: CM0007070



Introduction

Rabbit monoclonal antibody against cytochrome c oxidase subunit 1 (COX1). Validated for Western blotting, IHC-P, IF/ICC, and ELISA in human, mouse, and rat samples. Essential for studying mitochondrial respiratory chain function.

[Learn More](#)

Parameter	Value
Product Name	COX1 Rabbit mAb
Remarks/Alias	CoxI; COX1
Species	Mouse
Gene ID (Mouse)	17708
Immunogen	Synthetic peptide corresponding to amino acids 414-514 of human COX1 (NP_904330.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00397
Protein Weight	57 kDa
Shipping	Ice bag