

Anti-Cox15 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q7Kzn9

Item Number: CM0014705

Introduction

Rabbit polyclonal antibody targeting human COX15, validated for WB, IHC-P, IF/ICC, and ELISA. Ideal for heme A biosynthesis, mitochondrial function, and disease research.

[Learn More](#)

Parameter	Value
Product Name	COX15 Rabbit pAb
Remarks/Alias	MC4DN6; CEMCOX2; COX15
Species	Human
GeneID (human)	1355
GeneID	1355
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 150-250 of human COX15 (NP_510870.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7KZN9
Protein Weight	46kDa
Shipping	Ice bag