

Anti-Cox1/Ptgs1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P23219

Item Number: CM0029922



Introduction

Rabbit polyclonal antibody against COX1/PTGS1, validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Suitable for inflammation and prostaglandin biosynthesis studies.

[Learn More](#)

Parameter	Value
Product Name	COX1/PTGS1 Rabbit pAb
Remarks/Alias	COX1; COX3; PHS1; PCOX1; PES-1; PGHS1; PTGHS; PGG/HS; PGHS-1; COX1/PTGS1
Species	Human
GeneID (Human)	5742
GeneID	5742
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23219
Protein Weight	69kDa
Shipping	Ice bag