

Anti-Cox1/Ptgs1 Rabbit Monoclonal Antibody For Wb, Elisa - P23219

Item Number: CM0005166



Introduction

COX1/PTGS1 Rabbit Monoclonal Antibody validated for Western blot and ELISA, with human and mouse cross-reactivity. Ideal for studying prostaglandin biosynthesis, inflammatory response, and platelet activation.

[Learn More](#)

Parameter	Value
Product Name	COX1/PTGS1 Rabbit mAb
Remarks/Alias	COX1; COX3; PHS1; PCOX1; PES-1; PGHS1; PTGHS; PGG/HS; PGHS-1; COX1/PTGS1
Species	Human
Gene ID (Human)	5742
Gene ID	5742
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P23219
Protein Weight	69 kDa
Shipping	Ice bag