


Anti-Cox2/Ptgs2 Polyclonal Antibody For Wb, Ihc-P, Elisa - P35354

Item Number: CM0017522



Introduction

Rabbit polyclonal antibody against human COX2/PTGS2 (aa 375-604). Validated for WB, IHC-P, ELISA; cross-reacts with human, mouse, rat. Ideal for inflammation, prostaglandin biosynthesis, and related research.

[Learn More](#)

Parameter	Value
Product Name	Anti-COX2/PTGS2 Polyclonal Antibody
Remarks/Alias	COX2; COX-2; PHS-2; PGG/HS; PGHS-2; hCox-2; GRIPGHS; COX2/PTGS2
Species	Human
GeneID (Human)	5743
GeneID	5743
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 375-604 of COX2/PTGS2 (NP_000954.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35354
Protein Weight	69kDa
Shipping	Ice bag