


Anti-Cox2/Ptgs2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P35354

Item Number: CM0009050



Introduction

Rabbit monoclonal antibody targeting human COX2/PTGS2, validated for WB, IHC-P, and ELISA. Detects human, mouse, and rat COX2. Ideal for studying inflammation, prostaglandin synthesis, and pro-resolving mediators.

[Learn More](#)

Parameter	Value
Product Name	COX2/PTGS2 Rabbit mAb
Remarks/Alias	COX2; COX-2; PHS-2; PGG/HS; PGHS-2; hCox-2; GRIPGHS; COX2/PTGS2
Species	Human
GeneID (Human)	5743
GeneID	5743
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 505-604 of human COX2/PTGS2 (P35354).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35354
Protein Weight	69kDa
Shipping	Ice bag