


Anti-Cox2/Ptgs2 Rabbit Monoclonal Antibody For Wb - P35354

Item Number: CM0005173



Introduction

Rabbit monoclonal antibody against human COX2/PTGS2. Validated for WB, IP, IF-P, IHC-P, and ELISA. Cross-reacts with mouse and rat. Suitable for inflammation, cancer, and prostaglandin research.

[Learn More](#)

Parameter	Value
Product Name	COX2/PTGS2 Rabbit mAb
Remarks/Alias	COX2; COX-2; PHS-2; PGG/HS; PGHS-2; hCox-2; GRIPGHS
Species	Human
Gene ID (Human)	5743
Gene ID	5743
Immunogen	Synthetic peptide Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35354
Protein Weight	69kDa
Shipping	Ice bag