


Anti-Cox2/Ptgs2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q05769

Item Number: CM0024999



Introduction

Rabbit polyclonal antibody targeting mouse COX2/PTGS2. Validated for WB, IF/ICC, and ELISA. Suitable for studies on prostaglandin biosynthesis, inflammation, and resolution pathways. Immunogen: synthetic peptide within amino acids 101-210.

[Learn More](#)

Parameter	Value
Product Name	COX2/PTGS2 Rabbit pAb
Remarks/Alias	COX2; Cox-2; PES-2; PHS-2; Pghs2; TIS10; PGHS-2; PHS II; gripghs
Species	Mouse
GeneID (human)	5743
GeneID	19225
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 101-210 of mouse COX2/PTGS2 (NP_035328.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse
SWISS	Q05769
Protein Weight	69 kDa
Shipping	Ice bag