

Anti-Cox5A Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P20674

Item Number: CM0022173



Introduction

Rabbit polyclonal antibody targeting human COX5A, validated for knockout models. Suitable for WB, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Detects a 16 kDa mitochondrial protein. Ideal for oxidative phosphorylation research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] COX5A Rabbit pAb
Remarks/Alias	COX5A; COX; COX-VA; VA; cytochrome c oxidase subunit 5A; 5A
Species	Human
Gene ID (Human)	9377
Gene ID	9377
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human COX5A (NP_004246.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20674
Protein Weight	16kDa
Shipping	Ice bag