

Anti-Cox5A Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P20674

Item Number: CM0027234



Introduction

Rabbit polyclonal COX5A antibody validated in WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Suitable for mitochondrial complex IV research. SWISS: P20674.

[Learn More](#)

Parameter	Value
Product Name	COX5A Rabbit pAb
Remarks/Alias	VA; COX; COX-VA; MC4DN20; COX5A
Species	Human
GeneID (Human)	9377
GeneID	9377
Immunogen	Recombinant protein; recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human COX5A (NP_004246.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20674
Protein Weight	17kDa
Shipping	Ice bag