


Anti-Cox5B Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P10606

Item Number: CM0005860



Introduction

Mouse monoclonal antibody targeting human COX5B (Cytochrome c oxidase subunit 5B, mitochondrial). Validated for WB, IF/ICC, IHC-P, ELISA. Cross-reacts with mouse and rat. Recombinant fusion protein immunogen. For oxidative phosphorylation research.

[Learn More](#)

Parameter	Value
Product Name	COX5B Mouse mAb
Remarks/Alias	COXVB
Species	Human
GeneID (Human)	1329
GeneID	1329
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-129 of human COX5B (NP_001853.2)
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10606
Protein Weight	14 kDa
Shipping	Ice bag