

Anti-Cox5B Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P10606

Item Number: CM0025308



Introduction

COX5B Rabbit pAb targets mitochondrial cytochrome c oxidase subunit 5B, validated for WB, IHC-P, IF/ICC, IP, ELISA. Reacts with human, mouse, rat. 14kDa. Ideal for oxidative phosphorylation research.

[Learn More](#)

Parameter	Value
Product Name	COX5B Rabbit pAb
Remarks/Alias	COXVB; COX5B
Species	Human
GeneID (Human)	1329
GeneID	1329
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 32-129 of human COX5B (NP_001853.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10606
Protein Weight	14kDa
Shipping	Ice bag