


Anti-Cox6A1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P12074

Item Number: CM0028157



Introduction

High-quality rabbit polyclonal antibody targeting COX6A1, a mitochondrial cytochrome c oxidase subunit. Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Supports oxidative phosphorylation research.

[Learn More](#)

Parameter	Value
Product Name	COX6A1 Rabbit pAb
Remarks/Alias	COX6A; CMTRID; COX6AL; COX6A1
Species	Human
Gene ID (Human)	1337
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12074
Protein Weight	12kDa
Shipping	Ice bag