

Anti-Cox6B1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P14854

Item Number: CM0004848



Introduction

Human COX6B1 rabbit monoclonal antibody validated for WB, IHC-P, and ELISA. Ideal for studying mitochondrial cytochrome c oxidase, oxidative phosphorylation, and mitochondrial diseases. Monoclonal, human-reactive, ~10kDa protein weight.

[Learn More](#)

Parameter	Value
Product Name	COX6B1 Rabbit mAb
Remarks/Alias	COXG; COX6B; MC4DN7; COXVIb1; COX6B1
Species	Human
GeneID (Human)	1340
GeneID	1340
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-86 of human COX6B1 (P14854).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P14854
Protein Weight	10kDa
Shipping	Ice bag