


## Anti-Cox6B1 Rabbit Polyclonal Antibody For Wb, Elisa - P14854

Item Number: CM0027702



### Introduction

Rabbit polyclonal antibody targeting human COX6B1 (cytochrome c oxidase subunit 6B1). Validated for WB and ELISA; reacts with human and mouse. Useful for mitochondrial oxidative phosphorylation studies.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | COX6B1 Rabbit pAb                                                                                                |
| Remarks/Alias    | COXG; COX6B; MC4DN7; COXVIb1; COX6B1                                                                             |
| Species          | Human                                                                                                            |
| GeneID (Human)   | 1340                                                                                                             |
| GeneID           | 1340                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-86 of human COX6B1 (NP_001854.1) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, ELISA                                                                                                        |
| Cross Reactivity | Human, Mouse                                                                                                     |
| SWISS            | P14854                                                                                                           |
| Protein Weight   | 10kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |