


Anti-Cox6B1 Rabbit Polyclonal Antibody For Wb, Elisa - P14854

Item Number: CM0014735



Introduction

Anti-COX6B1 Rabbit Polyclonal Antibody validated for WB and ELISA. Detects human COX6B1, cross-reacts with mouse and rat. Ideal for mitochondrial respiratory chain complex IV studies.

[Learn More](#)

Parameter	Value
Product Name	COX6B1 Rabbit pAb
Remarks/Alias	COXG; COX6B; MC4DN7; COXVIb1; COX6B1
Species	Human
Gene ID (Human)	1340
Gene ID	1340
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-86 of human COX6B1 (NP_001854.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P14854
Protein Weight	10kDa
Shipping	Ice bag