


Anti-Cox6C Polyclonal Antibody For Wb, If/Icc, Elisa - P09669

Item Number: CM0020015



Introduction

COX6C rabbit polyclonal antibody validated for Western blot, IF/ICC, and ELISA. Recognizes human and mouse COX6C, a mitochondrial cytochrome c oxidase subunit essential for oxidative phosphorylation.

[Learn More](#)

Parameter	Value
Product Name	COX6C Rabbit pAb
Remarks/Alias	COX6C
Species	Human
GeneID (Human)	1345
GeneID	1345
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-75 of human COX6C (NP_004365.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse
SWISS	P09669
Protein Weight	9kDa
Shipping	Ice bag