

Anti-Cox7A1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P24310

Item Number: CM0017237

Introduction

Rabbit polyclonal antibody against COX7A1 (cytochrome c oxidase subunit 7A1), suitable for WB, IHC-P, ELISA. Detects human, mouse, rat COX7A1. Ideal for studying mitochondrial respiratory complex IV and oxidative phosphorylation.

[Learn More](#)

Parameter	Value
Product Name	COX7A1 Rabbit pAb
Remarks/Alias	COX7A; COX7AH; COX7AM; COX7A1
Species	Human
Gene ID	1346 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-79 of human COX7A1 (NP_001855.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24310
Protein Weight	9kDa
Shipping	Ice bag