

Anti-Cox7A2 Rabbit Polyclonal Antibody For Wb And Elisa - P14406

Item Number: CM0026732



Introduction

Rabbit polyclonal antibody targeting human COX7A2 (cytochrome c oxidase subunit 7A2, COX VIIa-L), validated for WB and ELISA. Ideal for studying mitochondrial oxidative phosphorylation and complex IV assembly.

[Learn More](#)

Parameter	Value
Product Name	COX7A2 Rabbit pAb
Remarks/Alias	VIIAL; COX7AL; COX7AL1; COXVIIAL; COXVIIa-L; COX7A2
Species	Human
GeneID (human)	1347
GeneID	1347
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 33-115 of human COX7A2 (NP_001856.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P14406
Protein Weight	9kDa
Shipping	Ice bag