

Anti-Cox7A2L Monoclonal Antibody For Wb, Ihc-P, Elisa - O14548

Item Number: CM0010025

Introduction

KD-validated rabbit monoclonal antibody targeting human COX7A2L/SCAF1. Cross-reacts with mouse and rat. Suitable for western blot, IHC-P, and ELISA. Ideal for studying mitochondrial respirasome assembly and bioenergetics.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] COX7A2L Rabbit mAb
Remarks/Alias	EB1; SCAF1; SCAFI; SIG81; COX7AR; COX7RP
Species	Human
GeneID (Human)	9167
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14548
Protein Weight	13kDa
Shipping	Ice bag