

Anti-Cox7A2L Rabbit Polyclonal Antibody (Kd Validated) For Wb And Elisa - O14548

Item Number: CM0022357

Introduction

Rabbit polyclonal antibody targeting COX7A2L (SCAF1), validated for Western blot (WB) and ELISA. Involved in mitochondrial respiratory supercomplex assembly; suitable for human sample research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] COX7A2L Rabbit pAb
Remarks/Alias	EB1; SCAF1; SCAFI; SIG81; COX7AR; COX7RP; [KD Validated] COX7A2L
Species	Human
Gene ID (Human)	9167
Gene ID	9167
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-114 of human COX7A2L (NP_004709.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O14548
Protein Weight	13kDa
Shipping	Ice bag