

Anti-Cpa3 Polyclonal Antibody For Wb, Elisa - P15088

Item Number: CM0020202

Introduction

Polyclonal antibody against human CPA3 (Mast cell carboxypeptidase A3), raised in rabbit. Validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Suitable for mast cell biology and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	CPA3 Rabbit pAb
Remarks/Alias	MC-CPA; CPA3
Species	Human
GeneID (Human)	1359
GeneID	1359
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 110-417 of human CPA3 (NP_001861.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P15088
Protein Weight	49kDa
Shipping	Ice bag