

Anti-Cpa6 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8N4T0

Item Number: CM0016167

Introduction

Anti-CPA6 rabbit polyclonal antibody validated for WB, IHC-P, ELISA. Detects human, mouse, rat CPA6 (51kDa), a secreted carboxypeptidase involved in neuropeptide processing and linked to epilepsy. SWISS: Q8N4T0.

[Learn More](#)

Parameter	Value
Product Name	CPA6 Rabbit pAb
Remarks/Alias	CPAH; ETL5; FEB11; CPA6
Species	Human
Gene ID	57094
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 130-437 of human CPA6 (NP_065094.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N4T0
Protein Weight	51kDa
Shipping	Ice bag