

Anti-Cpb1 Rabbit Polyclonal Antibody For Wb, Elisa - P15086

Item Number: CM0014525



Introduction

Rabbit polyclonal antibody against human CPB1, validated for WB and ELISA, cross-reacting with mouse and rat. Suitable for pancreatic carboxypeptidase B research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	CPB1 Rabbit pAb
Remarks/Alias	CPB; PASP; PCPB; CPB1
Species	Human
Gene ID (Human)	1360
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 111-417 of human CPB1 (NP_001862.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15086
Protein Weight	47kDa
Shipping	Ice bag