

Anti-Cplx2 Polyclonal Antibody For Wb, Elisa - Q6Puv4

Item Number: CM0024842

Introduction

Rabbit polyclonal antibody against human CPLX2 (Complexin-2), validated for WB and ELISA. Reacts with human, mouse, rat. Suitable for synaptic vesicle exocytosis and neurological disease research.

[Learn More](#)

Parameter	Value
Product Name	CPLX2 Rabbit pAb
Remarks/Alias	CPX2; Hfb1; 921-L; CPX-2; CPLX2
Species	Human
GeneID (Human)	10814
GeneID	10814
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-134 of human CPLX2 (NP_006641.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6PUV4
Protein Weight	15kDa
Shipping	Ice bag