

Anti-Cpne3 Polyclonal Antibody For Wb And Elisa - 075131

Item Number: CM0017596

Introduction

Anti-CPNE3 Rabbit Polyclonal Antibody validated for Western blotting (WB) and ELISA. Recognizes human and mouse CPNE3, a calcium-dependent phospholipid-binding protein involved in ERBB2-mediated tumor cell migration. Ideal for cancer and cell signaling studies.

[Learn More](#)

Parameter	Value
Product Name	CPNE3 Rabbit pAb
Remarks/Alias	CPN3, PRO1071, CPNE3
Species	Human
GeneID (Human)	8895
GeneID	8895
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 238-537 of human CPNE3 (NP_003900.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	075131
Protein Weight	60kDa
Shipping	Ice bag