

Anti-Cpox Rabbit Polyclonal Antibody For Wb And Elisa - P36551

Item Number: CM0028456



Introduction

Rabbit polyclonal antibody targeting human CPOX, validated for WB and ELISA applications. Cross-reacts with mouse and rat. Suitable for heme biosynthesis and mitochondrial research. Immunogen: recombinant CPOX (aa111-454).

[Learn More](#)

Parameter	Value
Product Name	CPOX Rabbit pAb
Remarks/Alias	COX; CPO; CPX; HCP; HARPO; CPOX
Species	Human
Gene ID	1371
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 111-454 of human CPOX (NP_000088.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36551
Protein Weight	50kDa
Shipping	Ice bag