


## Anti-Cps1 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P31327

Item Number: CM0010246



### Introduction

Rabbit monoclonal antibody against human CPS1, validated for WB, IF-P, IHC-P, ELISA. Targets mitochondrial urea cycle enzyme; removes excess ammonia. Ideal for liver and small intestine research.

[Learn More](#)

| Parameter        | Value                                                                                                                       |
|------------------|-----------------------------------------------------------------------------------------------------------------------------|
| Product Name     | CPS1 Rabbit mAb                                                                                                             |
| Remarks/Alias    | PHN; GATD6; CPSASE1; CPS1                                                                                                   |
| Species          | Human                                                                                                                       |
| Gene ID (Human)  | 1373                                                                                                                        |
| Gene ID          | 1373                                                                                                                        |
| Immunogen        | Synthetic Peptide   A synthetic peptide corresponding to a sequence within amino acids 100-200 of human CPS1 (NP_001866.2). |
| Source           | Rabbit                                                                                                                      |
| Category         | Monoclonal Antibodies                                                                                                       |
| Application      | WB, IF-P, IHC-P, ELISA                                                                                                      |
| Cross Reactivity | Human                                                                                                                       |
| SWISS            | P31327                                                                                                                      |
| Protein Weight   | 165kDa                                                                                                                      |
| Shipping         | Ice bag                                                                                                                     |