


Anti-Cps1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P31327

Item Number: CM0006612



Introduction

CPS1 rabbit monoclonal antibody for Western blot, IHC-P, and ELISA. Detects human, mouse, rat CPS1. Ideal for urea cycle and ammonia metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	CPS1 Rabbit mAb
Remarks/Alias	PHN; GATD6; CPSASE1; CPS1
Species	Human
GeneID	1373
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1356-1482 of human CPS1 (NP_001116105.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31327
Protein Weight	165kDa
Shipping	Ice bag