


Anti-Cps1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P31327

Item Number: CM0026068



Introduction

Rabbit polyclonal antibody targeting human CPS1 (CPSase I), validated for WB, IF/ICC, and ELISA. Detects CPS1 in human, mouse, and rat samples, suitable for urea cycle and ammonia detoxification research.

[Learn More](#)

Parameter	Value
Product Name	CPS1 Rabbit pAb
Remarks/Alias	PHN; GATD6; CPSASE1; CPS1
Species	Human
Gene ID (Human)	1373
Gene ID	1373
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 750-970 of human CPS1 (NP_001866.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31327
Protein Weight	165kDa
Shipping	Ice bag