


Anti-Cps1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P31327

Item Number: CM0020448



Introduction

High-quality rabbit polyclonal antibody against human CPS1 (carbamoyl-phosphate synthetase I). Validated for WB, IF/ICC, ELISA. Cross-reacts with rat. Immunogen: recombinant protein (aa 820-930). Ideal for urea cycle studies.

[Learn More](#)

Parameter	Value
Product Name	CPS1 Rabbit pAb
Remarks/Alias	PHN; GATD6; CPSASE1; CPS1
Species	Human
GeneID (Human)	1373
GeneID	1373
Immunogen	`Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	P31327
Protein Weight	165kDa
Shipping	Ice bag