


Anti-Cpsf1 Rabbit Polyclonal Antibody For Wb, Elisa - Q10570

Item Number: CM0024336



Introduction

High-quality rabbit polyclonal antibody targeting human CPSF1 (CPSF160), a cleavage and polyadenylation specificity factor subunit. Validated for Western blot and ELISA, with cross-reactivity to human and mouse. Ideal for studying pre-mRNA 3'-end formation and gene expression regulation.

[Learn More](#)

Parameter	Value
Product Name	CPSF1 Rabbit pAb
Remarks/Alias	MYP27; CPSF160; P/cl.18; HSU37012; CPSF1
Species	Human
GeneID (Human)	29894
GeneID	29894
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 700-950 of human CPSF1 (NP_037423.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q10570
Protein Weight	161kDa
Shipping	Ice bag