

Anti-Cpsf2 Polyclonal Antibody For Wb, Elisa - Q9P2I0

Item Number: CM0025882



Introduction

Rabbit polyclonal antibody targeting human CPSF2, validated in WB and ELISA with cross-reactivity to mouse and rat. Suitable for studying pre-mRNA 3'-end processing and histone mRNA maturation.

[Learn More](#)

Parameter	Value
Product Name	CPSF2 Rabbit pAb
Remarks/Alias	CPSF100; CPSF2
Species	Human
GeneID (Human)	53981
GeneID	53981
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 450-700 of human CPSF2 (NP_059133.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P2I0
Protein Weight	88kDa
Shipping	Ice bag