

Anti-Cpsf4L Polyclonal Antibody For Wb, Elisa - A6Nmk7

Item Number: CM0019713

Introduction

Rabbit polyclonal antibody targeting mouse Cpsf4l, validated for Western blot and ELISA. Suitable for studying mRNA processing and gene regulation in mouse samples. Recognizes a 21 kDa protein.

[Learn More](#)

Parameter	Value
Product Name	Cpsf4l Rabbit pAb
Remarks/Alias	D11ErtD636e; 0610010C04Rik; 1500000C01Rik; Cpsf4l
Species	Mouse
Gene ID (Human)	642843
Gene ID	52670
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-228 of mouse Cpsf4l (NP_080958.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	A6NMK7
Protein Weight	21kDa
Shipping	Ice bag