

Anti-Cpsf6 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q16630

Item Number: CM0026427



Introduction

Rabbit polyclonal antibody targeting CPSF6 (CFIm68). Validated for WB, IHC-P, IF/ICC, IP, ELISA; detects human, mouse, rat CPSF6 (59 kDa). Key factor in mRNA 3'-end processing and alternative polyadenylation.

[Learn More](#)

Parameter	Value
Product Name	CPSF6 Rabbit pAb
Remarks/Alias	CFIM; CFIM68; CFIM72; HPBR11-4; HPBR11-7; CPSF6
Species	Human
GeneID (human)	11052
GeneID	11052
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-330 of human CPSF6 (NP_008938.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16630
Protein Weight	59kDa
Shipping	Ice bag