

Anti-Cpsf73 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q9Ukf6

Item Number: CM0004691



Introduction

High-quality rabbit monoclonal antibody targeting CPSF73. Validated for western blot, immunoprecipitation, and ELISA; cross-reacts with human, mouse, and rat. Useful for studying mRNA 3'-end processing and histone mRNA cleavage.

[Learn More](#)

Parameter	Value
Product Name	CPSF73 Rabbit mAb
Remarks/Alias	CPSF73; NEDMHS; CPSF-73; NEDMHSN
Species	Human
Gene ID (Human)	51692
Gene ID	51692
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 585-684 of human CPSF3 (Q9UKF6).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UKF6
Protein Weight	77kDa
Shipping	Ice bag