

Anti-Cpt1A Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Ip, Elisa - P50416

Item Number: CM0010549



Introduction

High-quality rabbit monoclonal Anti-CPT1A antibody validated for WB, IF/ICC, IF-P, IHC-P, IP, ELISA. Reacts with human, mouse, rat. Ideal for studying fatty acid oxidation, mitochondrial uptake, and liver metabolism.

[Learn More](#)

Parameter	Value
Product Name	CPT1A Rabbit mAb
Remarks/Alias	CPT1; CPT1-L; L-CPT1
Species	Human
GeneID (Human)	1374
GeneID	1374
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 520-720 of human CPT1A (NP_001867.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50416
Protein Weight	88kDa
Shipping	Ice bag