

Anti-Cpt1A Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P50416

Item Number: CM0007068



Introduction

KO validated CPT1A Rabbit Monoclonal Antibody for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Targets liver isoform of carnitine palmitoyltransferase 1A, essential for fatty acid oxidation. Ideal for metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CPT1A Rabbit mAb
Remarks/Alias	CPT1; CPT1-L; L-CPT1; CPT1A
Species	Human
Gene ID (Human)	1374
Gene ID	1374
Immunogen	Recombinant protein of human CPT1A.
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50416
Protein Weight	88 kDa
Shipping	Ice bag