

Anti-Cpt1A Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P50416

Item Number: CM0030263



Introduction

Rabbit polyclonal antibody targeting human CPT1A, validated for WB, IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat CPT1A. Ideal for mitochondrial fatty acid oxidation studies.

[Learn More](#)

Parameter	Value
Product Name	CPT1A Rabbit pAb
Remarks/Alias	CPT1; CPT1-L; L-CPT1; CPT1A
Species	Human
Gene ID (Human)	1374
Gene ID	1374
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 520-720 of human CPT1A (NP_001867.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50416
Protein Weight	88kDa
Shipping	Ice bag