

Anti-Cpt1B Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q92523

Item Number: CM0025481

Introduction

CPT1B Rabbit Polyclonal Antibody (UniProt Q92523) validated for WB, IF-P, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for fatty acid metabolism and mitochondrial beta-oxidation research.

[Learn More](#)

Parameter	Value
Product Name	CPT1B Rabbit pAb
Remarks/Alias	CPT1; CPT1M; MCPT1; CPT1-M; CPTI-M; M-CPT1; MCCPT1; CPT1B
Species	Human
GeneID	1375
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 533-772 of human CPT1B (NP_004368.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92523
Protein Weight	88kDa
Shipping	Ice bag