

Anti-Cpt1B Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q92523

Item Number: CM0010124

Introduction

Anti-CPT1B rabbit monoclonal antibody validated for WB, IF-P, IHC-P, ELISA. Targets human, mouse, rat carnitine palmitoyltransferase 1B. Ideal for fatty acid metabolism and mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	CPT1B Rabbit mAb
Remarks/Alias	CPTI; CPT1M; MCPT1; CPT1-M; CPTI-M; M-CPT1; MCCPT1
Species	Human
Gene ID (Human)	1375
Gene ID	1375
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92523
Protein Weight	66 kDa/79-88 kDa
Shipping	Ice bag