

Anti-Cpt1C Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Tcg5

Item Number: CM0014169



Introduction

Rabbit polyclonal antibody against human CPT1C (Q8TCG5). Validated for WB, IF/ICC, ELISA with cross-reactivity to mouse and rat. Suitable for neuronal palmitoyl thioesterase and lipid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	CPT1C Rabbit pAb
Remarks/Alias	CATL1; CPT1P; CPTIC; SPG73; CPT1-B; CPTI-B; CPT1C
Species	Human
GeneID	126129
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 270-400 of human CPT1C (NP_001129524.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TCG5
Protein Weight	91kDa
Shipping	Ice bag