

Anti-Cpt2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P23786

Item Number: CM0017497

Introduction

Rabbit polyclonal antibody targeting human CPT2, validated for WB, IHC-P, IF/ICC, ELISA; cross-reacts with mouse, rat. Detects ~74kDa mitochondrial protein essential for fatty acid beta-oxidation.

[Learn More](#)

Parameter	Value
Product Name	CPT2 Rabbit pAb
Remarks/Alias	CPT1; IIAE4; CPTASE; CPT2
Species	Human
GeneID (Human)	1376
GeneID	1376
Immunogen	Recombinant protein corresponding to amino acids 404-658 of human CPT2 (NP_000089.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23786
Protein Weight	74kDa
Shipping	Ice bag