

Anti-Cpt2 Rabbit Monoclonal Antibody For Wb And Elisa - P23786

Item Number: CM0003273

Introduction

Anti-CPT2 rabbit monoclonal antibody validated for WB and ELISA in human, mouse, rat. Targets mitochondrial CPT2, essential for fatty acid beta-oxidation. Ideal for metabolism research.

[Learn More](#)

Parameter	Value
Product Name	CPT2 Rabbit mAb
Remarks/Alias	CPT1; IIAE4; CPTASE; CPT2
Species	Human
GeneID (Human)	1376
GeneID	1376
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 559-658 of human CPT2 (P23786).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23786
Protein Weight	74kDa
Shipping	Ice bag