

Anti-Cpt2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P23786

Item Number: CM0026765

Introduction

Rabbit polyclonal antibody against human CPT2, validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Suitable for fatty acid metabolism and mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	CPT2 Rabbit pAb
Remarks/Alias	CPT1; IIAE4; CPTASE; CPT2
Species	Human
GeneID	1376
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 420-658 of human CPT2 (NP_000089.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23786
Protein Weight	74kDa
Shipping	Ice bag