

Anti-Cpvl Rabbit Polyclonal Antibody For Wb, Elisa - Q9H3G5

Item Number: CM0013400

Introduction

Rabbit polyclonal antibody targeting human CPVL, validated for Western blot and ELISA. Detects mouse and rat orthologs. Ideal for studying lysosomal carboxypeptidase function in macrophages, phagocytosis, and antigen presentation.

[Learn More](#)

Parameter	Value
Product Name	CPVL Rabbit pAb
Remarks/Alias	HVLP; CPVL
Species	Human
GeneID (Human)	54504
GeneID	54504
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 127-476 of human CPVL (NP_061902.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9H3G5
Protein Weight	54kDa
Shipping	Ice bag