

Anti-Cpxm1 Polyclonal Antibody For Wb And Elisa - Q96Sm3

Item Number: CM0021908

Introduction

Rabbit polyclonal antibody against human CPXM1, validated for WB and ELISA. Detects human and mouse CPXM1, predicted molecular weight 82 kDa. Useful for studying cell-cell interactions and metalloprotease biology.

[Learn More](#)

Parameter	Value
Product Name	CPXM1 Rabbit pAb
Remarks/Alias	CPX1; CPXM; CPXM1
Species	Human
GeneID (Human)	56265
GeneID	56265
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 435-734 of human CPXM1 (Q96SM3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96SM3
Protein Weight	82kDa
Shipping	Ice bag