

Anti-Cr2/Cd21 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - P20023

Item Number: CM0003166



Introduction

High-quality rabbit monoclonal antibody against human CR2/CD21 (Complement receptor type 2). Validated for WB, IF-P, ELISA; reacts with human and rat. Immunogen: intracellular peptide (aa 934-1033). Ideal for B-cell signaling studies.

[Learn More](#)

Parameter	Value
Product Name	CR2/CD21 Rabbit mAb
Remarks/Alias	CR; C3DR; CD21; CVID7; SLEB9; CR2/CD21
Species	Human
GeneID (Human)	1380
GeneID	1380
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 934-1033 of human CR2/CD21 (P20023).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P20023
Protein Weight	113 kDa
Shipping	Ice bag